


IMPLANT SYSTEM


NEOBITECH CMI IMPLANT SYSTEM

- IS-II ACTIVE IMPLANT
- IS-III ACTIVE IMPLANT
- IS-III ACTIVE S-NARROW IMPLANT
- S-MINI IMPLANT
- SURGICAL KITS




FEATURES


IS-III Active Implant (S.L.A)



IS-II Active Implant (S.L.A)





Neobiotech Prosthetic Features




FEATURES

IS-III Active S-Narrow Implant (S.L.A)



S.L.A Surface

The new S.L.A. Surface with 40% greater fur-face area and 50% more cell adhesion promotes osseointegration



Deep & Wide Pitch

Optimum Pitch for Osseointegration (0.9Pitch)



Straight Body

Easy to adjust implantation depth



Wide Cutting Edge


Wide cutting edge and enlarged surface area enhances initial fixation and offers clinicians more stable implant placement




M2.0 Connection

- One Piece Type**
- Cover Screw
 - Healing Abutment

Internal Conical Seal 11° 2.1 Hex




Grade23 : Ti-6Al-4V ELI



M1.6 Connection Connection

- Two-Piece Type**
- PickCap Impression Coping
 - Temporary Abutment
 - Cement Abutment
 - Angled Abutment





IS-III Active Implant (S.L.A)

Length(mm)		7.3	8.5	10.0	11.5	13.0
Narrow	Ø3.2(P3.0)		IS3SN3008AP	IS3SN3010AP	IS3SN3011AP	IS3SN3011AP
Regular	Ø3.5(P3.7)		IS33508AP	IS33510AP	IS33511AP	IS33513AP
	Ø4.0(P4.1)	IS34007AP	IS34008AP	IS34010AP	IS34011AP	IS34013AP
	Ø4.5(P4.2)	IS34507AP	IS34508AP	IS34510AP	IS34511AP	IS34513AP
	Ø5.0(P4.35)	IS35007AP	IS35008AP	IS35010AP	IS35011AP	IS35013AP
S-Wide	Ø5.5(P4.35)	IS35507AP	IS35508AP	IS35510AP	IS35511AP	IS35513AP
	Ø6.0(P4.4)	IS36007AP	IS36008AP	IS36010AP	IS36011AP	IS36013AP

* Cover screw is included separately

Cover Screw


Type	Narrow	Regular / Wide / S-Wide
Prouct Code	ISN30S	ISC40S



IS-II Active Implant (S.L.A)

Length(mm)		7.3	8.5	10.0	11.5	13.0
Regular	Ø3.5(P3.5)		BIS3508A	BIS3510A	BIS3511A	BIS3513A
	Ø4.0(P3.7)	BIS4007A	BIS4008A	BIS4010A	BIS4011A	BIS4013A
	Ø4.5(P3.8)	BIS4507A	BIS4508A	BIS4510A	BIS4511A	BIS4513A
	Ø5.0(P4.0)	BIS5007A	BIS5008A	BIS5010A	BIS5011A	BIS5013A
S-Wide	Ø5.5(P4.53)	BIS5507A	BIS5508A	BIS5510A	BIS5511A	BIS5513A
	Ø6.0(P4.9)	BIS6007A	BIS6008A	BIS6010A	BIS6011A	BIS6013A
	Ø7.0(P5.8)	BIS7007A	BIS7008A	BIS7010A	BIS7011A	BIS7013A

* Cover screw is included with Fixture (Narrow/All 7.3mm Fixtures: ISN30S, Regular/Wide/S-Wide Fixtures: ISC40S)




Healing Abutment

- Hand Tightening

Cuff(mm)		2.0	3.0	4.0	5.0	6.0	7.0
Diameter	Ø4.0	ISUH4002	ISUH4003	ISUH4004	ISUH4005	ISUH4006	
	Ø4.5	ISUH4502	ISUH4503	ISUH4504	ISUH4505	ISUH4506	
	Ø4.8	ISH402	ISH403	ISH404	ISH405	ISH406	ISH407
	Ø5.5	ISH502	ISH503	ISH504	ISH505	ISH506	ISH507
	Ø6.0	ISH602	ISH603	ISH604	ISH605	ISH606	ISH607
	Ø6.8	ISH702	ISH703	ISH704	ISH705	ISH706	ISH707
	Ø8.0	ISH802	ISH803	ISH804			
	Ø9.0	ISH902	ISH903	ISH904			

* Ø8.0, Ø9.0 Healing Abutments are for Custom Abutment




ABUTMENT

- Abutment Screw(ISCS20) included, Tightening Torque: 30Ncm


IS Cemented Abutment (Hex Type)

Cuff(mm)			1.0	2.0	3.0	4.0	5.0
Diameter	Ø4.5	4.0mm	ISAH4140	ISAH4240	ISAH4340	ISAH4440	ISAH4540
		5.5mm	ISAH415	ISAH425	ISAH435	ISAH445	ISAH455
		7.0mm	ISAH417	ISAH427	ISAH437	ISAH447	ISAH457
	Ø5.2	4.0mm	ISAH5140	ISAH5240	ISAH5340	ISAH5440	ISAH5540
		5.5mm	ISAH515	ISAH525	ISAH535	ISAH545	ISAH555
		7.0mm	ISAH517	ISAH527	ISAH537	ISAH547	ISAH557
	Ø5.7	4.0mm	ISAH6140	ISAH6240	ISAH6340	ISAH6440	ISAH6540
		5.5mm	ISAH615	ISAH625	ISAH635	ISAH645	ISAH655
		7.0mm	ISAH617	ISAH627	ISAH637	ISAH647	ISAH657
	Ø6.5	4.0mm	ISAH7140	ISAH7240	ISAH7340	ISAH7440	ISAH7540
		5.5mm	ISAH715	ISAH725	ISAH735	ISAH745	ISAH755




IS Cemented Abutment (SCRIP Type)

Cuff(mm)			1.0	2.0	3.0	4.0	5.0
Diameter	Ø4.5	4.0mm	ISAS4140	ISAS4240	ISAS4340	ISAS4440	ISAS4540
		5.5mm	ISAS415	ISAS425	ISAS435	ISAS445	ISAS455
		7.0mm	ISAS417	ISAS427	ISAS437	ISAS447	ISAS457
	Ø5.2	4.0mm	ISAS5140	ISAS5240	ISAS5340	ISAS5440	ISAS5540
		5.5mm	ISAS515	ISAS525	ISAS535	ISAS545	ISAS555
		7.0mm	ISAS517	ISAS527	ISAS537	ISAS547	ISAS557
	Ø5.7	4.0mm	ISAS6140	ISAS6240	ISAS6340	ISAS6440	ISAS6540
		5.5mm	ISAS615	ISAS625	ISAS635	ISAS645	ISAS655
		7.0mm	ISAS617	ISAS627	ISAS637	ISAS647	ISAS657
	Ø6.5	4.0mm	ISAS7140	ISAS7240	ISAS7340	ISAS7440	ISAS7540
		5.5mm	ISAS715	ISAS725	ISAS735	ISAS745	ISAS755



IS Cemented Abutment (Non-Hex Type)

Cuff(mm)			1.0	2.0	3.0	4.0	5.0
Diameter	Ø4.5	4.0mm	ISAN4140	ISAN4240	ISAN4340	ISAN4440	ISAN4540
		5.5mm	ISAN415	ISAN425	ISAN435	ISAN445	ISAN455
		7.0mm	ISAN417	ISAN427	ISAN437	ISAN447	ISAN457
	Ø5.2	4.0mm	ISAN5140	ISAN5240	ISAN5340	ISAN5440	ISAN5540
		5.5mm	ISAN515	ISAN525	ISAN535	ISAN545	ISAN555
		7.0mm	ISAN517	ISAN527	ISAN537	ISAN547	ISAN557
	Ø5.7	4.0mm	ISAN6140	ISAN6240	ISAN6340	ISAN6440	ISAN6540
		5.5mm	ISAN615	ISAN625	ISAN635	ISAN645	ISAN655
		7.0mm	ISAN617	ISAN627	ISAN637	ISAN647	ISAN657
	Ø6.5	4.0mm	ISAN7140	ISAN7240	ISAN7340	ISAN7440	ISAN7540
		5.5mm	ISAN715	ISAN725	ISAN735	ISAN745	ISAN755



ABUTMENT

• Abutment Screw(ISC20) included, Tightening Torque: 30Ncm


IS Angled Abutment

15 degree

Cuff(mm)		2.0	3.0	4.0
Ø4.5	A Type	ISAHA1427	ISAHA1437	ISAHA1447
	B Type	ISAHB1427	ISAHB1437	ISAHB1447
	Non-Hex	ISANA1427	ISANA1437	ISANA1447
Ø5.2	A Type	ISAHA1527	ISAHA1537	ISAHA1547
	B Type	ISAHB1527	ISAHB1537	ISAHB1547
	Non-Hex	ISANA1527	ISANA1537	ISANA1547
Ø5.7	A Type	ISAHA1627	ISAHA1637	ISAHA1647
	B Type	ISAHB1627	ISAHB1637	ISAHB1647
	Non-Hex	ISANA1627	ISANA1637	ISANA1647

25 degree

Cuff(mm)		2.0	3.0	4.0
Ø4.5	A Type	ISAHA2427	ISAHA2437	ISAHA2447
	B Type	ISAHB2427	ISAHB2437	ISAHB2447
	Non-Hex	ISANA2427	ISANA2437	ISANA2447
Ø5.2	A Type	ISAHA2527	ISAHA2537	ISAHA2547
	B Type	ISAHB2527	ISAHB2537	ISAHB2547
	Non-Hex	ISANA2527	ISANA2537	ISANA2547
Ø5.7	A Type	ISAHA2627	ISAHA2637	ISAHA2647
	B Type	ISAHB2627	ISAHB2637	ISAHB2647
	Non-Hex	ISANA2627	ISANA2637	ISANA2647




* Customizable by grinding. Can be positioned in 12 directions by selecting A or B type

UCLA Abutment

Gold UCLA


Hex	Non-Hex
ISUGH400	ISUGN400



UCLA Abutment

Plastic UCLA

Diameter	Hex	Non-Hex
Ø4.5	ISUPH400	ISUPN400
Ø5.2	ISUPH500	ISUPN500
Ø5.7	ISUPH600	ISUPN600



TruScan Body

Type	Platform Size	Length(mm)	Product Code
IS	Regular	One Size	DUR-SB
	Mini (S-Narrow only)		OTM-SB
IT	Regular	One Size	STR-SB
	Wide (Ø6.5 & above)		STW-SB




Tru Base

Type	Platform Size	Length(mm)	Product Code
IS	Regular	3.5	DUR-B13S
		5.5	DUR-B15S
	Mini (S-Narrow only)	3.5	OTM-B13S
		5.5	OTM-B15S
IT	Regular Neck	3.5	STR-B13S
	Wide Neck (Ø6.5 & above)		STW-B13S



IS Temporary Abutment

Cuff(mm)	2.0	2.0	2.0
Length(mm)	6.0	8.0	13.0
Hex Type Ø4.5	ISAHT460	ISAHT480	ISAHT4130
SCRP Type Ø4.5	ISAST460	ISAST480	ISAST4130
Non-Hex Type Ø4.5	ISANT480	ISANT480	ISANT4130




MULTI UNIT ABUTMENT - NEOBIOTECH

• Screw-retained abutment for multiple or bridges, Tightening Torque: 30Ncm, Cylinder Screw: 15Ncm

IS Multi Unit Abutment


				
Cuff(mm)	1.0	2.0	3.0	4.0
Product Code	ISMUA4810S	ISMUA4820S	ISMUA4830S	ISMUA4840S



* Holder(MUAH48) is included

IS Multi Unit Angled Abutment


					
Cuff(mm)		1.0	2.0	3.0	4.0
Hex	17°	ISAHMUA1410S	ISAHMUA1420S	ISAHMUA1430S	ISAHMUA1440S
	30°	ISAHMUA3410S	ISAHMUA3420S	ISAHMUA3430S	ISAHMUA3440S
Non-Hex	17°	ISANMUA1410S	ISANMUA1420S	ISANMUA1430S	ISANMUA1440S
	30°	ISANMUA3410S	ISANMUA3420S	ISANMUA3430S	ISANMUA3440S



* Holder is included

Temporary Cylinder, Non-Hex


Length(mm)	11.0
Product Code	MUATPCN48



* Cylinder Screw(MUAS 14) is included

Transfer Impression Coping

Length(mm)	9.0
Product Code	MUAIT48




Multi Unit Abutment Driver

Product Code	MUAD48S
--------------	---------




Healing Cap

Length(mm)	4.1
Product Code	MUAHC48



Pick-up Impression Coping

Length(mm)	9.0
Product Code	MUAIP48




Lab Analog

Product Code	MUALA48
--------------	---------



CCM Cylinder, Non-Hex

Length(mm)	12.0
Product Code	MUACCN48



* Cylinder Screw(MUAS 14) is included

1.2 Hex Driver

		
Length(mm)	10.0	15.0
Product Code	HD1210S	HD1215S

KERATOR ABUTMENT

Regular Fixture (Ø3.5~Ø7.0)

Straight Abutment

Set		Straight Abutment						Male Cap				
Cuff(mm)		1.0	2.0	3.0	4.0	5.0	6.0					
Diameter	Ø4.0	AS401	AS402	AS403	AS404	AS405	AS406	Metal Housing	White Spacer	2.4 lbs	1.2 lbs	Extra Light Angled



Angled Abutment

Set	Angled Abutment		Male Cap				
Diameter	Ø4.0						
Degree	15°	15°					
Cuff(mm)	1.5	3.0					
Product Code	AS40AN1.5	AS40AN3.0					



S-Narrow Fixture (Ø3.2)

Straight Abutment

Set		Straight Abutment						Male Cap				
Product Code		G351	G352	G353	G354	G355	G356					
								Metal Housing	White Spacer	2.4 lbs	1.2 lbs	Extra Light Angled

Male Processing Cap

				
Metal Housing	White Spacer	2.4 lbs	1.2 lbs	Extra Light Angled
Product Code		CDPHSET		

Male Cap (Each)

Type (Color)	Product Code
White	CDPH WHITE
Black	CDPH BLACK
Blue	CDPH BLUE
Red	CDPH RED
Pink	CDPH PINK



Tool

Product Code	KMT002
--------------	--------



Tip

Product Code	KMD719
--------------	--------



Coping

Product Code	DKI4845
--------------	---------



Analog

Product Code	DKA3854
--------------	---------




SOLID ABUTMENT & RELATED

- Abutment Screw(ISC520) included, Tightening Torque: 30Ncm

IS Solid Abutment

Cuff(mm)		1.0	2.0	3.0	4.0	5.0
Ø4.5	4.0mm	ISAE4140	ISAE4240	ISAE4340	ISAE4440	ISAE4540
	5.5mm	ISAE415	ISAE425	ISAE435	ISAE445	ISAE455
	7.0mm	ISAE417	ISAE427	ISAE437	ISAE447	ISAE457
Ø5.2	4.0mm	ISAE5140	ISAE5240	ISAE5340	ISAE5440	ISAE5540
	5.5mm	ISAE515	ISAE525	ISAE535	ISAE545	ISAE555
	7.0mm	ISAE517	ISAE527	ISAE537	ISAE547	ISAE557
Ø5.7	4.0mm	ISAE6140	ISAE6240	ISAE6340	ISAE6440	ISAE6540
	5.5mm	ISAE615	ISAE625	ISAE635	ISAE645	ISAE655
	7.0mm	ISAE617	ISAE627	ISAE637	ISAE647	ISAE657
Ø6.5	5.5mm	ISAE715	ISAE725	ISAE735	ISAE745	ISAE755



Protective Cap

Diameter Ø4.5

Length(mm)	4.5	5.5	7.0
Product Code	ISPTC445	ISPTC455	ISPTC470

Diameter Ø5.7


Length(mm)	4.5	5.5	7.0
Product Code	ISPTC645	ISPTC655	ISPTC670

Diameter Ø5.2

Length(mm)	4.5	5.5	7.0
Product Code	ISPTC545	ISPTC555	ISPTC570

Diameter Ø6.5

Length(mm)	4.5	5.5	7.0
Product Code	ISPTC745	ISPTC755	ISPTC770



Impression Cap

Diameter Ø4.5

Product Code	ISPIC411
--------------	----------

Diameter Ø5.7


Product Code	ISPIC611
--------------	----------

Diameter Ø5.2

Product Code	ISPIC511
--------------	----------


Diameter Ø6.5

Product Code	ISPIC711
--------------	----------



Lab Analog

Length(mm)	4.0	5.5	7.0
Ø4.5	ISLA440	ISLA455	ISLA470
Ø5.2	ISLA540	ISLA555	ISLA570
Ø5.7	ISLA640	ISLA655	ISLA670
Ø6.5	ISLA740	ISLA755	ISLA770



Plastic Coping

Diameter	Ø4.5	Ø5.2	Ø5.7	Ø6.5
Hex	ISPCH410	ISPCH510	ISPCH610	ISPCH710
Non-Hex	ISPCN410	ISPCN510	ISPCN610	ISPCN710



BALL ABUTMENT

- Abutment Screw(ISCS20) included, Tightening Torque: 30Ncm
- Compensation degree: Max 20° between implants

Ball Abutment


Cuff(mm)	1.0	2.0	3.0	4.0
Product Code	ISABA100	ISABA200	ISABA300	ISABA400



Housing Set

Product Code	BAH40
--------------	-------


- O-ring(BAORING) 5ea. and O-ring(BAOIMP) 1ea. Included



Retainer Set


Product Code	BAR20
--------------	-------

- O-ring(BAORING) 5ea. and O-ring(BAOIMP) 1ea. Included



O-Ring (Clinical)


Product Code	BAORING
--------------	---------



O-Ring (Lab)

Product Code	BAOIMP
--------------	--------

- 3ea. Included
- Prevents the resin flowing into the under of housing or Retainer
 - While making a final resin denture at laboratory
 - While inserting Housing or Retainer within denture directly at the office




Ball Lab Analog

Product Code	BALA350
--------------	---------



Ball Abutment Driver


Product Code	BADH24
--------------	--------



IMPRESSION COPING

Pick-up Impression Coping (Open Tray)

Diameter	Length(mm)	Hex	SCRP	Non-Hex
Ø4.0	Short(12.0mm)	ISIPH400	ISIPS400	ISIPN400
	Long(15.0mm)	ISIPHL400	ISIPSL400	ISIPNL400
Ø4.5	Short(12.0mm)	ISIPH450	ISIPS450	ISIPN450
	Long(15.0mm)	ISIPHL450	ISIPSL450	ISIPNL450
Ø4.8	Short(12.0mm)	ISIPH411	ISIPS411	ISIPN411
	Long(15mm)	ISIPHL411	ISIPSL411	ISIPNL411
Ø5.5	Short(12.0mm)	ISIPH511	ISIPS511	ISIPN511
	Long(15.0mm)	ISIPHL511	ISIPSL511	ISIPNL511
Ø6.0	Short(12.0mm)	ISIPH611	ISIPS611	ISIPN611
	Long(15.0mm)	ISIPHL611	ISIPSL611	ISIPNL611




Impression Coping Screw




Short	ISPIS20
Long	ISPISL20

Transfer Impression Coping (Closed Tray)

Diameter	Length(mm)	Hex	Non-Hex
Ø4.8	Short(12.0mm)	ISITH411	ISITN411
	Long(15.0mm)	ISITHL411	ISITNL411
Ø5.5	Short(12.0mm)	ISITH511	ISITN511
	Long(15.0mm)	ISITHL511	ISITNL511
Ø6.0	Short(12.0mm)	ISITH611	ISITN611
	Long(15.0mm)	ISITHL611	ISITNL611




Impression Coping HEX Screw



Short	ISTIS20H
Long	ISTISL20H

Lab Analog



Diameter	Ø3.2	Ø3.8	Ø4.3
Product Code	ISLA300	ISLA400	ISLA500


Clinical Guide

Fixture	Ø3.2	Ø3.5	Ø4.0 and above
Analog	ISLA300 Only	ISLA400 Only	ISLA500 Only
Abutment	Ø3.5 Only	All abutment except Ø3.5	All abutment except Ø3.5

* S-Narrow implants (Ø3.2mm) refer to page18.

IS-III ACTIVE S-NARROW FIXTURE

Fixture



Diameter	Length(mm)	Product Code
Ø3.2	8.5	IS3SN3008AP
Ø3.2	10.0	IS3SN3010AP
Ø3.2	11.5	IS3SN3011AP
Ø3.2	13.0	IS3SN3013AP


* Cover Screw is included separately

Cover Screw



Diameter	Product Code
Ø3.0	ISC25S


Healing Abutment



Diameter	Ø3.0	Ø3.0
Cuff(mm)	3.0	5.0
Product Code	ISH3003	ISH3005


Diameter	Ø3.8	Ø3.8
Cuff(mm)	3.0	5.0
Product Code	ISH3803	ISH3805

Temporary Abutment




Diameter	Length(mm)	Connection	Product Code
Ø3.2	12.0	Hex	ISAHT300
		Non-Hex	ISANT300

Cement Abutment (Hex)




Diameter	Post(mm)	Cuff(mm)	Product Code
Ø3.5	7.0	3.0	ISAH3537
		4.0	ISAH3547
		5.0	ISAH3557
		7.0	ISAH3577

Cement Abutment (Non-Hex)



Diameter	Post(mm)	Cuff(mm)	Product Code
Ø3.5	7.0	3.0	ISAN3537
		4.0	ISAN3547
		5.0	ISAN3557
		7.0	ISAN3577

PickCap Impression Coping



Diameter	Length(mm)	Product Code
Ø3.5	5.0	ISSNPCI3505C
	7.0	ISSNPCI3507C

* Screw & Cap are included

Cap Size	Product Code	Remark
Small	PCISC0110S(10ea.)	Below Ø4.5

Lab Analog




Diameter	Product Code
Ø3.2	ISLA300

S-MINI FIXTURE

S-MINI CEMENT TYPE (RBM SURFACE)


Length(mm)	8.5	10.0	11.5	13.0	15.0
Ø2.5	MI2508L	MI2510L	MI2511L	MI2513L	MI2515L
Ø3.0	MI3008L	MI3010L	MI3011L	MI3013L	MI3015L
Ø3.5	MI3508L	MI3510L	MI3511L	MI3513L	MI3515L

* Protective Cap is Included



S-MINI BALL TYPE (RBM SURFACE)

Length(mm)	8.5	10.0	11.5	13.0	15.0
Ø3.0	MIB30083	MIB30103	MIB30113	MIB30133	MIB30153
Ø3.5	MIB35083	MIB350103	MIB35113	MIB35133	MIB35153



Impression Cap

Product Code	MICIMP
--------------	--------




Protective Cap

Product Code	MICPTC
--------------	--------



Lab Analog - Cement Type

Type	Normal	Modifying
Product Code	MICLAB01	MICLAB02



Lab Analog - Ball Type


Product Code	BALA350
--------------	---------



Housing Set

Product Code	BAH40
--------------	-------


- O-ring(BAORING) 5ea. and O-ring(BAOIMP) 1ea. Included



O-ring (Clinical)

Product Code	BAORING
--------------	---------


- 3ea. Included



Retainer Set

Product Code	BAR20
--------------	-------


- In case of occlusion interval is lower, use retainer instead of housing
- O-ring(BAORING) 5ea. and O-ring(BAOIMP) 1ea. Included




O-ring (Lab)

Product Code	BAOIMP
--------------	--------

- 3ea. Included
- Prevents the resin flowing into the under of Housing or Retainer
 - while making a final resin denture at laboratory.
 - while inserting Housing or Retainer within denture directly at the office.



NAVI GUIDE KIT




Navi Guide Kit

Product Code	NEONGKIT
--------------	----------


Tissue Punch

Diameter(Ø)	Product Code
Ø4.3	NGTP43R




Point Drill

Diameter(Ø)	Product Code
Ø2.2	NGPD22




Bone Trimmer

Diameter(Ø)	Type	Product Code
Ø2.2	Standard	NGBT00
	1.5mm Offset	NGBT15
	3.0mm Offset	NGBT30




Ø2.2 Initial Taper Drill

Diameter(Ø)	Type	Product Code
Ø2.2	7.3	NGTWD2207
	8.5	NGTWD2208
	10.0	NGTWD2210
	11.5	NGTWD2211
	13.0	NGTWD2213




Stop Drill




Diameter(Ø)		Ø3.0	Ø3.5	Ø4.0	Ø4.5
Length (mm)	7.3	NGTD3007	NGTD3507	NGTD4007	NGTD4507
	8.5	NGTD3008	NGTD3508	NGTD4008	NGTD4508
	10.0	NGTD3010	NGTD3510	NGTD4010	NGTD4510
	11.5	NGTD3011	NGTD3511	NGTD4011	NGTD4511
	13.0	NGTD3013	NGTD3513	NGTD4013	NGTD4513

NEO PICKCAP IMPRESSION KIT



Product Code	PCIKIT
--------------	--------

Neo PickCap Impression Mini Kit Tray




Product Code	PCIKITTB / KTC
--------------	----------------




Neo PickCap Impression Coping

Cuff(mm)		2.0	3.0	4.0	5.0	6.0	8.0	10.0
Diameter	Ø3.5		ISPCIB3503C		ISPCIB3505C		ISPCIB3508C	ISPCIB3510C
	Ø4.0		ISPCIB4003C		ISPCIB4005C		ISPCIB4008C	ISPCIB4010C
	Ø4.5		ISPCIB403C		ISPCIB405C		ISPCIB408C	ISPCIB410C
	Ø5.2	ISPCIB502C	ISPCIB503C	ISPCIB504C	ISPCIB505C	ISPCIB506C		
	Ø5.7	ISPCIB602C	ISPCIB603C	ISPCIB604C	ISPCIB605C	ISPCIB606C		
	Ø6.5	ISPCIB702C	ISPCIB703C	ISPCIB704C	ISPCIB705C	ISPCIB706C		
	Ø8.0		ISPCIB803C	ISPCIB804C				
	Ø9.0		ISPCIB903C	ISPCIB904C				



Neo PickCap Impression Coping Cap

Cap Size	Product Code	Remark
Small	PCISC0110S(10ea.)	Below Ø4.5
Large	PCILC0110S(10ea.)	Over Ø5.2




Neo PickCap Impression Coping Screw

Product Code	ISPCIS20
--------------	----------



Hex Driver

Length(mm)	Product Code
15.0	HD1215S
20.0	HD1220S
25.0	HD1225S



PickCap Positioner

Length(mm)	Product Code
2.0	AP202N





Pincette



Product Code	PCIPC
--------------	-------


STOP DRILL KIT



Stop Drill

Drill Diameter		Ø2.2	Ø3.0	Ø3.5	Ø4.0	Ø4.5	Ø5.0
Length (mm)	7.3	TSD2207S9	TSD3007S9	TSD3507S9	TSD4007S9	TSD4507S9	TSD5007S9
	8.5	TSD2208S9	TSD3008S9	TSD3508S9	TSD4008S9	TSD4508S9	TSD5008S9
	10.0	TSD2210S9	TSD3010S9	TSD3510S9	TSD4010S9	TSD4510S9	TSD5010S9
	11.5	TSD2211S9	TSD3011S9	TSD3511S9	TSD4011S9	TSD4511S9	TSD5011S9
	13.0	TSD2213S7	TSD3013S7	TSD3513S7	TSD4013S7	TSD4513S7	TSD5013S7

* OPTIONAL: Drill Ø5.0(7.3~13.0)



Guide Point Drill


Diameter	Product Code
Ø1.35	PD13



Lindermann Drill

Diameter	Product Code
Ø2.3	LDS23C

* OPTIONAL, Lindermann Drill is NOT included in Stop Drill Kit



STOP DRILL KIT

Cortical Tap (0.8(IS-II)/0.9(IS-III) Pitch)

		Cortical Tap (IS-II Active)		Cortical Tap (IS-III Active)	
Diameter	FØ3.5	CT38S08		CT38S09	
	FØ4.0	CT43S08		CT43S09	
	FØ4.5	CT45S08		CT45S09	
	FØ5.0	CT50S08		CT50S09	
	FØ5.5	CT55S08		CT55S09	

* OPTIONAL: Cortical Tap Ø5.5/IS-II Active, Ø5.5/IS-III Active

IS Fixture Driver

Length	Product Code
Ratchet(Short)	ISFD10R
Ratchet(Long)	ISFD15R
Contra Angle(Short)	ISFD05C
Contra Angle(Long)	ISFD05CL




MD Parallel Pin

Diameter	Product Code
Ø4.0 / Ø5.0	MDPP4050
Ø6.0 / Ø7.0	MDPP6070
Ø8.0 / Ø9.0	MDPP8090



Direction Pin

Diameter	Product Code
Ø3.6	DP36
Ø4.5	DP45




Torque Ratchet

Product Code	TW60
--------------	------



Connector

Length	Product Code
Long	RC15
Product Code	CAA00



Hex Driver

Length(mm)	Product Code
15.0	HD1215S



IT Fixture Driver

Length	Product Code
Ratchet(Short)	ITFD031RS
Ratchet(Long)	ITFD031RL
Contra Angle(Short)	ITFD031S
Contra Angle(Long)	ITFD031L



* OPTIONAL, IT Fixture Driver is NOT included in Stop Drill Kit

MD Spacer

Diameter	Product Code
Ø4.0	MDS4025
Ø7.0	MDS7030
Ø8.0	MDS8030



* Only for Guide Point Drill

Drill Extension

Product Code	DE01
--------------	------




Driver Holder & Depth Gauge

Product Code	DHDG
--------------	------




IS CORE KIT




Product Code	NEOCKIT
--------------	---------

Drilling Sequence (IS-III Active)



Drilling Sequence (IS-II Active)



IS CORE KIT

Guide Point Drill

Diameter	Product Code
Ø2.3	PD13




Surgical Drill

Diameter	Product Code
Ø2.2	TSD22CS
Ø3.0 (3.5F)	TSD30CS
Ø3.5 (4.0F)	TSD35CS
Ø4.0 (4.5F)	TSD40CS
Ø4.5 (5.0F)	TSD45CS



Cortical Tap (IS-III Active Implant)

Diameter	Product Code
Ø3.5	CTS38S09
Ø4.0	CTS43S09
Ø4.5	CTS45S09
Ø5.0	CTS50S09



Drill Stopper

					
Length(mm)	7.3	8.5	10.0	11.5	13.0
Product Code	DS070C	DS085C	DS100C	DS115C	DS130C


IS Fixture Driver

Product Code	ISFD05C
--------------	---------



Parallel Pin

Length(mm)	Product Code
8.5	PP08F
10.0	PP10F




Drill Extension

Product Code	DE01
--------------	------



Ratchet Connector

Product Code	RC15
--------------	------



Hex Driver

Product Code	HD1215S
--------------	---------




Torque Ratchet

Product Code	TW60
--------------	------




S-WIDE KIT




S-WIDE TWIST SURGICAL Drill

Diameter	Product Code
Ø4.8	TSD48
Ø5.0	TSD50
Ø5.5	TSD55
Ø6.0	TSD60
Ø6.5	TSD65




S-WIDE Multi Profile

Diameter	Product Code
Ø5.5	MP55
Ø6.0	MP60
Ø7.0	MP70




Trephine Drill

Diameter	Product Code
Ø4.0 / Ø5.0	TB40
Ø5.0 / Ø6.0	TB50



S-WIDE Surgical Tap

Diameter	Product Code
Ø5.5	TD55
Ø6.0	TD60
Ø7.0	TD70




S-WIDE Drill Stopper

Length(mm)	6.0	7.0	8.5	10.0	11.5
Product Code	DSL060	DSL070	DSL085	DSL100	DSL115


HEX Driver

Type	Product Code
1.2 Hex	HD1212s



Connector

Type	Product Code
Ratchet	RC10
Contra Angle	CAA00



S-NARROW KIT



Component	Product Code
KIT TRAY TOP	CSTKITTRAY1
KIT TRAY MIDDLE	SNKITTRAY2
KIT TRAY BOTTOM	CSTKITTRAY3

* The top, middle, bottom tray and components of the S-narrow Kit are sold separately

Drill Extension



Product Code
DE01

Cortical Tap



Diameter	Product Code
Ø3.2	CT32S09

Hex Driver



Product Code
HD1215S

Guide Point Drill



Diameter	Product Code
Ø1.35	PD13

Cortical Drill



Diameter	Product Code
Ø3.1	ISCD32

Direction Pin



Diameter	Product Code
Ø3.6	DP36

Surgical Drill



Diameter	Product Code
Ø2.2	TSDS22CL
Ø2.7	TSD27CL

Fixture Driver



Product Code
ISFDSNRL
ISFDSNCL

Connector



Product Code
RC15S

Drill Stopper




Length(mm)	8.5	10.0	11.5	13.0
Product Code	DS085CS	DS100CS	DS115CS	DS130CS

Torque Ratchet




Product Code	TW60
--------------	------

S-MINI KIT




S-MINI Guide Drill

Product Code	MGD12
--------------	-------



S-Mini Twist Surgical Drill


Diameter	Product Code
Ø1.5	MSD15
Ø2.0	MSD20
Ø2.5	MSD25
Ø3.0	MSD30



S-Mini Surgical Tap

Diameter	Product Code
Ø2.5	MTD25**
Ø3.0	MTD30**
Ø3.5	MTD35**


** OPTIONAL



Connector

Length	Product Code
Ratchet(Short)	RC10**
Ratchet(Long)	RC15
Contra Angle	CAA00

** OPTIONAL



Fixture Driver

Type	Product Code
Cement	ISQD20
Ball	BADH24




Parallel Pin

Product Code	MPP10
--------------	-------



IMPLANT SYSTEM



NEOBIOTECH USA

2700 E. FOOTHILL BLVD. #121
PASADENA, CA 91107
www.neobiotechusa.com
213-387-7701, 7704